

TO REACH
PROFESSIONALS IN LABORATORY SCIENTISTS
WORLDWIDE

ADVERTISE
ONLINE



LABPLUSINTERNATIONAL.COM

The online source of technology & product
information for analytical & laboratory scientists

[Magazine](#) | [This Week's News](#) | [E-products](#) | [Events](#) | [Contact us](#) | [My Lab-on](#)



Cell culture CO2 incubator

Cell culture CO2 incubator



With a unique H₂O₂ decontamination option that limits downtime to less than three hours, the new MCO-19AIC CO₂ incubator is ideal when total chamber decontamination with verification is required. The H₂O₂ decontamination process functions via the company's SafeCell UV system. Following a seven minute H₂O₂ vaporisation, circulation and dwell cycle, the vaporisation is stopped and the UV lamp is turned on for 90 minutes. The new single beam, dual detector infrared CO₂ sensing system offers unprecedented control accuracy and stability by simultaneously measuring sample and reference wavelengths for continuous auto-zero calibration..

Sanyo E&E Europe BV

Nijverheidsweg 120
Po Box 494 - 4870 AL
Etten-Leur 4879 AZ
NETH
Tel.: +31- 76-5433838
Fax: +31- 76-5413732



Lab
International

ADVERTISE
ONLINE



CLICK
HERE
TO FIND
OUT MORE

